

Corning Incorporated
Life Sciences

Registered ISO 9001

Catalog Number CLIP-01
PDS reference PDS_Cbor_2262 Version _ 2
Product description Corning® Gosselin™ Clip for Blender Bag, Aseptic, 50/Bag, 200/Box, 1000/Case

Sanitary Assurance Level

ASEPTIC products are manufactured in a controlled environment, without human intervention during the manufacturing process. Finished products with an Aseptic claim are not subjected to sterilization and do not make a claim to have a Sterility Assurance Level.

Expiration Statement

N/A

Performance testing (AQL)

Inspection based on a statistically significant number of samples towards to internal or agreed specifications that determine Assurance Quality Level.

Primary body composition

PE

Secondary body composition

N/a

Closure

N/a

Special feature

Aseptic

Packaging

50/Bag, 200/Box, 1000/Case



Hugues Vanderschelden
Plant Quality Manager, Borre
vandersch2@corning.com