

Catalog Number: 3732

Product Description: Corning® CellSTACK® - One Piece Universal Cap

**Component Materials:**

- One Piece Universal Cap - Virgin High-Density Polyethylene, meets USP, Class VI requirements for plastic Containers and closures. Heavy metal free (meets CONEG req.) color concentrate.
- Over cap - Low Density Polyethylene
- Vent - 0.2µm PTFE membrane, meets USP, Class VI requirements for plastic containers and closures.
- Foil - Polyethylene heat sealed film

**Product Dimensions:**

- |                 |   |              |                |   |          |
|-----------------|---|--------------|----------------|---|----------|
| Diameter of Cap | - | 1.46 in.     | Height of Cap  | - | 1.5 in.  |
| Tolerances      | - | +/- 0.05 in. | Overall Height | - | 1.60 in. |

**Sterilization** – Product has been sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137 “Sterilization of health care products- Radiation”. Products meet a minimum Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

**Quality Control Testing** – Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release. Key inspections and inline tests are listed below:

- Visual Inspection – Pass
- Packaging Inspection – Pass
- Leak Test – Pass