



## ROTOR SPECIFICATIONS

Type	Fixed Angle
Material	Titanium
Capacity	20 x 0.5 mL
Maximum Speed	70,000 rpm
Tube Dimensions Ø x L	7 x 56 mm
Angle	30°
K Factor	30.7
Net Weight	1.4 kg/ 3.1 lbs
Warranty	5 years

	RCF	RADIUS
Maximum	307,052 xg	56.1 mm
Average	238,636 xg	43.6 mm
Minimum	169,672 xg	31.0 mm

## ROTOR AND ACCESSORIES ORDERING INFORMATION

## S70-AT Rotor Package includes:

Description	Cat. No.
S70-AT Rotor	45591
0.5 mL PC Thick-Walled Tubes (qty. 1 pkg. of 100)	45243
Large O-ring for Cover (qty. 2)	45297
Small O-ring for Cover (qty. 2)	45294
Rotor Holding Fixture	45266
Lubricant	61556
Vacuum Grease	65937
Instruction Manual	45735

## Accessories:

Description	Cat. No.
Rotor Stand	45357
Cover Opener	45173

## LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
0.5	0.5	PC Thick-Walled Tube	45243	100	20	70,000	307,052	7 x 56	—	—	—	—

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit  
[www.thermoscientific.com/rotors](http://www.thermoscientific.com/rotors).

DSCFGMULTRAS70AT 1114

**Thermo**  
 SCIENTIFIC

Part of Thermo Fisher Scientific