



The following features are present in the vector based on the nucleotide sequence.

Multiple cloning site (MCS)	1-42
Red Firefly luciferase gene	49-1695
BGH poly(A) signal	1707-1832
SV40 origin/promoter	1833-2397
Puromycin resistant gene	2398-2997
SV40 poly(A) signal	3159-3192
beta-lactamase (Amp ^R) gene	3301-4161
pUC replication origin (pUC Ori)	4340-5144
Transcriptional terminator (Ter)	5145-5752
Lac operator 1 (Lac O1)	5753-5773
Transcriptional pause site (TPS)	5906-5977