

## **TECHNICAL SHEET**

PRODUCT NAME: Anti-Giardia mAb cyst protein (clone GR16)

Code: MT-16GR16

Clone: GR16

Isotype: IgG<sub>1</sub>

Physical State: Liquid

Host Animal (source): Mouse

**Specificity:** GR16 is specific for Giardia detection.

Cross Reactivity: No cross reactivity with Cryptosporidium parvum, Entamoeba

histolytica, Entamoeba dispar.

**Purification:** Anti-Giardia mAb cyst protein (clone GR16) has been purified in a

multi-step procedure. ≥95% pure by SDS-PAGE.

**Storage Conditions:** Storage at 2-8°C.

**Health & Safety Information:** This product contains sodium azide as preservative. Although the

amount is 0.09%, this substance is poisonous and hazardous and

it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance

with the respective national regulations.

**Presentation:** PBS pH: 7.4, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 11/01/2019

