

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 8502

**Product Description:** Corning® 96 1D/2D Bar Coded Storage Tube, 1.3 mL, With Screw Cap, Polypropylene, Not Sterile, Bulk Packed in Bag, 96 Tubes/Bag, 10 Bags/Case

**Component Materials:**

Tube - Polypropylene  
Cap - Polypropylene

**Attributes:**

Working Volume - 1000 µL  
Product Dimensions - 43.9 Outside Length w/o Cap; 42.5 mm Inside Depth; 6.8 mm Opening ID  
Cap Dimensions - 8.7 x 14.3 mm

**Heating / Cooling:**

This product is designed for use at temperatures ranging from -100°C (148°F) to 110°C (230°F).

**Special Features:**

Each tube features a unique identification number in 1D bar code format (Code 128) and human readable text on the side of the tube and in 2D bar code format (14x14 Data Matrix Code with ECC200 error correction) on the bottom of the tube.

**Sterilization:**

Supplied Not Sterile.

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date; Next 2 digits - Last 2 digits of year of manufacture; Last 3 digits - Batch identification.

Rev No: 2