

## Ordering Information

Cat#		CKS-A103Z050	CKS-A105T050	CKS-A106Z080	CKS-A108Z080	CKS-A110T125	CKS-A115T216	CKS-A120Q216	CKS-A125Q216	
Unit parameter										
Geometric capacity (L)		3	5	7.5	8	10	15.9	21	25	
Neck diameter(mm)		50	50	80	80	125	216	216	216	
Static evaporation loss rate (L/Day)		0.16	0.17	0.20	0.22	0.43	0.84	0.84	0.84	
Static holding time (Day)		20	31	37	35	23	18	25	29	
Storage capacity										
Canisters	Canister outer diameter (mm)	38	38	63	63	97	—	—	—	
	Canister height (mm)	120	276	120	120	276	—	—	—	
	Canisters (ea)	1	1	1	1	1	—	—	—	
	Straw holding capacity (single layer)	0.5ml(ea)	132	214	820	820	1508	—	—	—
		0.25ml(ea)	298	472	1780	1780	3324	—	—	—
Racks and vial boxes	Racks (ea)	—	—	—	—	1	1	1	1	
	Box size (ea)	—	—	—	—	76×76	134×134	134×134	134×134	
	Layers per racks (ea)	—	—	—	—	4	3	4	5	
	Vial (ea)	—	—	—	—	100	300	400	500	
25ml	Racks (ea)	—	—	—	—	1	1	1	1	
	Layers of racks (ea)	—	—	—	—	2	1	2	3	

