

Thermo Scientific SuperSealer

Ensure a proper seal on plates, blocks and storage tube racks with the Thermo Scientific SuperSealer automatic sealer. Choose individual sealing options such as Thermo Scientific Matrix® SepraSeal® mats for liquid handling and single tube accessibility or seal entire plates and blocks for storage with Thermo Scientific Matrix CapMats. Effectively apply adhesive tape seals and ensure a tight seal with the included SuperSealer Adhesive Seal Mat.



Sealing Simplicity. The Thermo Scientific SuperSealer simplifies the sealing of plates, blocks and tube racks of virtually any height, with no adjustments needed, ensuring a secure and consistent seal. Once the seal along with the plate, block or tube rack are placed into the drawer, simply close the door to initiate sealing.

Sealing Flexibility. Designed to work with Thermo Scientific Matrix SepraSeal mats, CapMats and, with the SuperSealer Adhesive Seal Mat, any pre-cut adhesive tape seal, you can be assured a proper seal regardless of sealing option.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate the SuperSealer in your laboratory, detailing how it can address your labware sealing needs.

Thermo Scientific SuperSealer

Speed	10-15 seconds/plate, block or rack of tubes*
Dimensions (W x D x H)	6.5 "/16.5 cm x 12.5 "/31.75 cm x 13.5 "/34.3 cm
Weight	30 lbs/13.5 kg
Power Supply	110-240 V

*The SuperSealer is designed to only work with full racks of tubes

Ordering Guide

Item No.	Description	Qty
4110	SuperSealer	1 each
4111	SuperSealer Adhesive Seal Mat	1 each

Compatible Thermo Scientific Matrix Sealing Options

Item No.	Description	Qty
4410	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks, Bulk	100/case
4411	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks	10/case
4412	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks, Sterile	10/case
4420	CapMat for 2 ml Blocks, Bulk	100/case
4421	CapMat for 2 ml Blocks	10/case
4422	CapMat for 2 ml Blocks, Sterile	10/case
4430	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks, Bulk	100/case
4431	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks	10/case
4432	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks, Sterile	10/case
4463	SepraSeal Mats	10 mats of 96 caps/case
4464	SepraSeal Mats, Sterile	10 mats of 96 caps/case
4465	SepraSeal Mats, Pre-Split	10 mats of 96 caps/case
4466	SepraSeal Mats, Sterile, Pre-Split	10 mats of 96 caps/case
4469	SepraSeal Cap Removal Tool	1 each
4417	Adhesive Sealing Tape for Plates and Blocks	100 seals/case
4418	Adhesive Sealing Tape for Plates and Blocks, DMSO-Resistant	100 seals/case
4419	Adhesive Sealing Tape for Plates and Blocks, Foil	100 seals/case