



**Comments for pYES7L:
7,266 nucleotides**

TRP1 gene: bases 172–846 [c]

ARS4/CEN5: bases 1,382–1,902

OriT: bases 2,196–2,307

rrnB T2: bases 2,450–2,477 [c]

rrnB T1: bases 2,609–2,652 [c]

trfA gene: bases 3,408–4,553 [c]

Spectinomycin resistance gene: bases 5,085–6,094 [c]

Spn (Spectinomycin promoter): bases 6,096–6,229 [c]

OriV: bases 6,243–6,859

[c] = complementary strand