

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 8505

Product Description: Corning® 96 Storage Tube Screw Cap With O-Ring, Polypropylene, Not Sterile, On Mat with Lid, 96 Caps/Mat, 5 Mats/Pack, 10 Packs/Case (4800 Caps)

Component Materials:

Cap - Polypropylene
Mat - Polypropylene
Lid - Polypropylene

Attributes:

Cap Dimensions - 8.7 mm OD x 14.3 mm Height
Rack (Mat) Dimensions - 82.0 x 118.7 mm Base; 21.6 mm Height (w/o Lid); 32.6 mm Height (w/Lid)
Hole Pattern - 8 x 12 on 9.0 mm centers

Heating / Cooling:

This product is designed for use at temperatures ranging from -100°C (148°F) to 110°C (230°F).

Sterilization:

Supplied Not Sterile.

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date; Next 2 digits - Last 2 digits of year of manufacture; Last 3 digits - Batch identification.

Rev No: 1