

Ordering Information

Unit parameter								
Cat#	CKS-P302	CKS-P302	CKS-P002	CKS-P603	CKS-P606	CKS-P608Z	CKS-P610RZ050	
	Z030	Z035	R050	Z050	Z050	080		
Geometric capacity (L)	2	2	2	3	6	8	10	
Neck diameter (mm)	30	35	50	50	50	80	50	
Static evaporation loss rate (L/D)	0.07	0.08	0.10	0.12	0.12	0.21	0.12	
Static holding time(Da)	28	24	20	26	52	38	86	
Storage capacity								
Canister outer diameter (mm)	19	25	—	38	38	63	38	
Number of canister (ea)	3	3	—	6	6	6	6	
Straw holding capacity (single layer)	0.5ml (ea)	90	165	—	792	792	2244	792
	0.25ml (ea)	204	330	—	1788	1788	5022	1788
Straw holding capacity (double layers)	0.5ml (ea)	—	—	—	—	—	—	—
	0.25ml (ea)	—	—	—	—	—	—	—
Containers								
Height (mm)	399	428	385	435	482	502	552	
Outer diameter (mm)	223	204	223	223	300	300	300	
Empty weight(kg)	2.8	2.6	2.7	3.1	4.8	5.9	5.9	
Accessories included								
Lockable lid	√	√	√	√	√	√	√	
Protective holding bag	√	√	√	√	√	√	√	
Optional accessories								
Intelligent monitor	√	√	√	√	√	√	√	