

Ordering Information

Cat#	CKST620Z/T050	CKST630Z/T050	CKST630Z/T080	CKST635Z/T050	CKST635Z/T080	CKST650Z/T050	CKST650Z/T080
Unit Parameter							
Geometric capacity (L)	20	31.5	31.5	35.5	35.5	50	50
Neck diameter(mm)	50	50	80	50	80	50	80
Static evaporation loss rate (L/Day)	0.20	0.20	0.30	0.20	0.30	0.23	0.34
Static holding time (Day)	101	159	106	179	119	213	147
Storage capacity							
Canister outer diameter (mm)	38	38	63	38	63	38	63
Canister height(mm)	120/276	120/276	120/276	120/276	120/276	120/276	120/276
Canisters (ea)	6	6	6	6	6	6	6
Straw holding capacity (single layer)	792	792	2244	792	2244	792	2244
	1788	1788	5022	1788	5022	1788	5022
Straw holding capacity (double layers)	1284	1284	3624	1284	3624	1284	3624
	2832	2832	8460	2832	8460	2832	8460
Containers							
Height (mm)	672	706	710	750	754	811	823
Outer diameter (mm)	394	462	462	462	462	462	462
Empty weight (kg)	10.1	13.3	13.4	14.4	14.8	15.4	15.7
Accessories included							
Lockable lid	√	√	√	√	√	√	√
Protective holding	√	√	√	√	√	√	√
Optional accessories							
Intelligent monitor	√	√	√	√	√	√	√
Carrying cart	—	√	√	√	√	√	√