


Corning® Cell Expansion Bag

CORNING


Corning's single-use Cell Expansion Bag features a single-web polyolefin film. This Cell Expansion Bag is intended for cell culture and is an ideal environment for the scalable expansion of suspension cell lines. The integrated tubing allows for closed system filling, feeding, and sampling. The recommended fill volume ranges from 190 mL to 1252 mL in four standard sizes.

Features

- ▶ Cell expansion observed with multiple cell models
- ▶ Gas permeable film
- ▶ Reusable sampling valve for in process testing
- ▶ Excellent film clarity when filled for culture viewing
- ▶ Tubing for sterile weld connections
- ▶ Scalable, user-friendly design
- ▶ Available in the following sizes: 500 mL, 1 L, 3 L, and 5 L




Ordering Information



Port Design - Thin walled membrane port, sampling valve, two PVC leads, female luers, male luer caps, and small pinch clamps.

Cat. No.	Description	Film Type	Unit Size	Fill Volume	Max Fill	Qty
91-200-84	Cell Expansion Bag	Polyolefin	500 mL	190 mL	510 mL	1
91-200-85	Cell Expansion Bag	Polyolefin	1 L	381 mL	1400 mL	1
91-200-86	Cell Expansion Bag	Polyolefin	3 L	633 mL	3280 mL	1
91-200-87	Cell Expansion Bag	Polyolefin	5L	1252 mL	5000 mL	1



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.cellgro.com or call 1.800.235.5476.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

Mediatech, Inc.
A Corning Subsidiary
 9345 Discovery Boulevard
 Manassas, VA 20109
 t 800.235.5476
 t 703.471.5955
 f 703.467.9851
www.cellgro.com

CORNING | **FALCON**® **AXYGEN**® **PYREX**® **GOSSSELIN**™

Corning cellgro products are manufactured by Mediatech, Inc.
 For a listing of Corning trademarks, visit us at www.cellgro.com/about-us/trademarks.html
 All other trademarks in this document are the property of their respective owners.