

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 26797-38

**Product Description:** PYREX® Closure, Polypropylene, Natural, 38 mm

**Component Materials:**  
Polypropylene (PP)

**Attributes:**  
Product Dimensions - 29 mm OD x 31 mm High

**Autoclaving:**  
This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

**Heating / Cooling:**  
Polypropylene may be heated to 140°C.

Polypropylene may be cooled to -196°C.

**Vacuum / Pressure:**  
This product was not designed for and is not rated for vacuum applications.

This product was not designed for and is not rated for pressure applications.

**Sterilization:**  
Shipped Not Sterile

**Lot Number Designation:**  
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**  
5 AUG 2013

Rev No: 1