

Catalog Number: 3567

Product Description: Corning® Disposable Spinner Flask Orange Vent Cap, GL45 thread, 0.2um membrane

Component Materials:

- Cap - Virgin High Density Polyethylene, meets *USP Class VI requirements for plastic containers and closures*. Heavy metal free (meets *CONEG* req.) color concentrate
- Filter - 0.2 µm microporous PTFE membrane, meets *USP Class VI requirements for plastic containers and closures*.

Product Dimensions:

Diameter	-	45mm (1.77 in.)	Height	-	25.3mm (1.0 in.)
Tolerances	-	+/- .050 in			

Animal Content:

Product does not contain materials of animal origin.

Sterility – Product has been sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137, “Sterilization of health care products- Radiation”. Products meet a minimum Sterility Assurance Level (SAL) of 10⁻⁶.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual and microscopic examination of the product.

Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits – Julian date, start of manufacturing; Next 2 digits – Year of manufacture; Last 3 digits – Batch identification.

Rev 5