

Compatibility table – pipette tips

Pipette tips with filter

2.5 µl		70.3010.265 70.3010.285	70.3010.365 70.3010.385
10 µl		70.3010.255 70.3010.275	70.3010.355 70.3010.375
20 µl		70.3020.255 70.3020.275	70.3020.355 70.3020.375
20 µl		70.3030.265 70.3030.285	70.3030.365 70.3030.385
100 µl		70.3030.255 70.3030.275	70.3030.355 70.3030.375
200 µl		70.1189.215	
200 µl		70.3031.255 70.3031.275	70.3031.355 70.3031.375
300 µl		70.3040.255 70.3040.275	70.3040.355 70.3040.375
1,000 µl		70.3050.255 70.3050.275	70.3050.355 70.3050.375
1,000 µl		70.3060.255 70.3060.275	70.3060.355 70.3060.375
1,250 µl		70.1186.210 70.1186.410	

Pipette tips without filter

10 µl		70.3010 70.3010.200 70.3010.210 70.3010.205	70.3010.305 70.3010.100 70.3010.110
20 µl		70.3020 70.3020.200 70.3020.205 70.3020.305	70.3020.210 70.3020.100 70.3020.110
20 µl		70.3021 70.3021.200 70.3021.205 70.3021.305	70.3021.210 70.3021.100 70.3021.110
200 µl		70.3030 70.3030.020 70.3030.200 70.3030.205	70.3030.305 70.3030.210 70.3030.100 70.3030.110
200 µl		70.1189.105	
250 µl		70.3031 70.3031.020 70.3031.200 70.3030.205	70.3031.305 70.3031.210 70.3031.100 70.3031.110
300 µl		70.3040 70.3040.100 70.3040.110 70.3040.200	70.3040.205 70.3040.210 70.3040.305
1,000 µl		70.3050 70.3050.020 70.3050.030 70.3050.100 70.3050.110	70.3050.200 70.3050.205 70.3050.210 70.3050.305
1,250 µl		70.3060 70.3060.100 70.3060.110 70.3060.200	70.3060.205 70.3060.210 70.3060.305
1,250 µl		70.1186 70.1186.100 70.1186.200	

Eppendorf® pipettes					Thermo® Finnpipettes					Gilson® pipettes																																		
0.1-2.5 µl	0.5-10 µl	2-20 µl	10-100 µl	50-200 µl	0.1-2.5 µl	0.5-10 µl	2-20 µl	10-100 µl	20-200 µl	0.1-2.5 µl	0.5-10 µl	2-20 µl	10-100 µl	20-200 µl	0.2-2 µl	F2 0.5-5 µl	F1 1-10 µl	0.5-10 µl Univ.	2-20 µl Univ.	5-50 µl MC	10-100 µl Univ.	20-200 µl Univ.	50-300 µl MC	100-1,000 µl	200-1,000 µl	0.3-3 µl	0.5-5 µl	1-10 µl	30-300 µl	30-300 µl	30-300 µl MC	100-1,000 µl	100-1,200 µl MC	50-300 µl MC	P20/20 0.2-2 µl	P100/100 0.1-10 µl	C8/C12 1-10 µl	P200/200 2-20 µl	P1000/1000 20-200 µl	P2000/2000 20-200 µl	U8 20-300 µl	P1000/1000 200-1,000 µl	P5000/5000 1,000-5,000 µl	P1000/1000 1-10 ml
Eppendorf Reference®	Eppendorf Research®	Eppendorf Research® plus	Eppendorf Research® pro	Xplore®	Finnpipette® Digital / F1 / F2	Finnpipette® Focus	Finnpipette® Novus	Bio Control	Pipetman® / Ultra / Concept																																			
● 2.5	● 2.5	● 2.5	● 2.5	● 2.5	● 2.5	● 2.5	● 2.5		● 2.5																																			
● 10	● 10	● 10	● 10	● 10	● 10	● 10	● 10		● 10																																			
● 20	● 20	● 20	● 20	● 20	● 20	● 20	● 20		● 20																																			
● 20	● 20	● 20	● 20	● 20	● 20	● 20	● 20		● 20																																			
● 100	● 100	● 100	● 100	● 100	● 100	● 100	● 100		● 100																																			
● 200	● 200	● 200	● 200	● 200	● 200	● 200	● 200		● 200																																			
● 200	● 200	● 200	● 200	● 200	● 200	● 200	● 200		● 200																																			
● 300	● 300	● 300	● 300	● 300	● 300	● 300	● 300		● 300																																			
● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000		● 1,000																																			
● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000	● 1,000		● 1,000																																			
● 1,250	● 1,250	● 1,250	● 1,250	● 1,250	● 1,250	● 1,250	● 1,250		● 1,250																																			

● = compatible
10 = reduced pipetting volume
X = not recommended
 = not checked
 MC = multi-channel pipette

Pipette manufacturers revise or modify the pipette cones on their pipettes as part of technical changes. Our recommendations in the compatibility table can therefore only be taken as a guide. In case of doubt, we are happy to provide sample tips to you free of charge in order to ensure the optimal fit and performance of our tips with your pipettes.

Technical modifications reserved
 This document may contain information on products that may not be available in particular countries
 20_809_0300_207

