CORNING

ITS (Insulin, Transferrin, and Selenium Solution) Formulation

Cat. No.	25-800-CR
Description	Liquid Solution, 100x
Units	mg/L
Components	
Selenious acid	0.67
Insulin, human recombinant	1000
Transferrin, human recombinant	550



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.cellgro.com or call 1.800.235.5476.

Mediatech, Inc. A Corning Subsidiary 9345 Discovery Boulevard Manassas, VA 20109 t 800.235.5476 t 703.471.5955 f 703.467.9851

www.cellgro.com

CORNING | FALCON° AXYGEN. **GOSSELIN**[™] **PYREX**°