

Lighting Incubator

Features:

- Microprocessor PID control, accurate and stable.
- Magnetic door sealing.
- LED display for BJPX-L series, LCD display for BJPX-LII series.
- CFC-free refrigerant, eco-friendly and highly efficient.

Safety Functions:

- Over temperature protection.
- Temperature sensor for open circuit or short circuit protection.
- Parameters can be saved automatically when power off or system error.



BJPX-L150/II BJPX-L250/II
BJPX-L300/II BJPX-L400/II

BJPX-L450II

Optional Accessories: Printer (only for BJPX-LII series)

Notice: Please don't keep setting the temperature lower than 20°C for more than 160h.

Technical Parameters:

Model	BJPX-L150/II	BJPX-L250/II	BJPX-L300/II	BJPX-L400/II	BJPX-L450II
Capacity	159L	248L	290L	402L	450L
Lighting Type	BJPX-L series: 2-side illumination; BJPX-LII series:3-side illumination				
Temp. Range	With illumination: 10-50°C, Without illumination: 5-50°C				
Temp. Fluctuation	±0.5°C				±1°C
Temp. Uniformity	±2°C				
Lighting Illumination	BJPX-L 150/250/300/400: 0-6000 Lux (5-grade adjustable);				
	BJPX-L150II/250II/300II/400II: 0-12000 Lux (5-grade adjustable);				
	BJPX-450II: 0-15000 (6-grade adjustable)				
Ambient Temp.	5-35°C				
Continuous Working Time	≥180h				
Refrigerant	R134a (CFC free, eco-friendly)				
Shelves	3 pos/adjustable				4 pos/adjustable
Consumption	1000W	1500W	2000W	2500W	1750W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)				
Internal Size(W*D*H) mm	420*420*928	470*449*1178	490*464*1278	560*600*1198	700*550*1200
External Size(W*D*H) mm	621*752*1575	672*792*1825	685*780*1920	755*905*1840	1030*780*1900
Package Size(W*D*H) mm	710*860*1740	770*890*1990	790*900*2100	870*1050*2040	1140*940*2080
Gross Weight(kg)	148	158	178	200	185