

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 6120P-50

**Product Description:** Corning® Funnel, Reusable Plastic, PP, Not Sterile, 50 mm

**Component Materials:**  
Polypropylene (PP)

**Attributes:**  
Product Dimensions - 50 mm Top ID, 11.5 - 8 mm Stem ID, 78 mm High

**Autoclaving:**  
Autoclaving this product using a standard cycle of 121°C and 15 psi for 15 minutes may change the physical shape of the product and its original volumetric accuracy.

**Heating / Cooling:**  
The upper temperature limit for this plastic material is 125°C.

This product IS NOT suitable for Dry Heat Sterilization as it will melt. This product MUST NOT be heated by any direct, surface contacting heat source.

A lower limit usage temperature has not been established for this product.

**Vacuum / Pressure:**  
This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**  
Shipped Not Sterile

**Lot Number Designation:**  
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**  
01 OCT 2012

Rev No: 1