

## High Temperature Drying Oven



### Features:

- \* Mirror stainless steel inner chamber.
- \* Optional independent temperature protector.
- \* Optional test hole diameter 50mm on the left side.
- \* Optional RS485 interface for software connection.

### Technical Parameters:

Model	BOV-H50F	BOV-H90F	BOV-H216F
Capacity	50L	100L	200L
Temp. Range	50~400°C		
Temp. Precision	0.1°C		
Temp. Fluctuation	±0.5°C		
Ambient Temp.	5~40°C		
Timing Range	1~9999min		
Shelves No.	2 pcs		
Consumption	2500W	3000W	3500W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)		380V, 50HZ
Internal Size(W*D*H) mm	355*380*400	450*450*480	600*600*610
External Size(W*D*H) mm	810*640*640	900*710*710	1060*860*860
Package Size(W*D*H) mm	940*800*800	1030*870*870	1190*1020*1020
Gross Weight(kg)	70	110	150