

Drying Oven / Incubator(Dual-use)



BOV-D50



BOV-D30



BOV-D70/BOV-D140/BOV-D240

Features:

- * Over-temperature protection.
- * Microprocessor temperature control with LED display.
- * Stainless steel inner chamber, round angle structure, adjustable shelves.
- * High quality aluminum silicate silk cotton six sides insulation
- * Leakage protector

Technical Parameters:

Model	BOV-D30	BOV-D50	BOV-D70	BOV-D140	BOV-D240
Capacity	30L	50L	70L	140L	240L
Temp. Range	Incubator: RT+5-80°C; Drying Oven: 80-200°C				
Temp. Precision	0.1°C				
Temp. Fluctuation	±0.5°C				
Ambient Temp.	5-40°C				
Timing Range	1-9999min				
Shelves No.	2 pcs				
Consumption	550W	850W	1100W	1500W	2050W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)				
Internal Size(W*D*H) mm	340*325*290	420*350*350	400*375*470	450*557*520	500*615*720
External Size(W*D*H) mm	490*520*690	710*520*520	540*550*870	600*740*920	650*800*1130
Package Size(W*D*H) mm	620*630*830	810*640*600	670*660*1040	730*860*1080	780*920*1270
Gross Weight(kg)	55	50	75	85	110