

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Product Description

Catalog Number: 6520

Product Description: Costar® 96-well Lid

Component Materials:

Plate - Polycarbonate

Product Dimensions:

Length of Lid	-	5.03 in.	Depth of Lid	-	0.50 in.
Width of Lid	-	3.40 in.	Tolerances	-	+/- 0.05 in.

Functional Characteristics:

The product is thermally stable up to 135°C, and can be used on Thermowell polycarbonate plates as additional protection to oil overlays.

Compatible PCR Plates:

Catalog numbers 6510, 6511, 6514.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1