



The following features are present in the vector based on the nucleotide sequence.

Cytomegalovirus (CMV) promoter.....	8-635
<i>Gaussia</i> luciferase gene	647-1204
BGH poly(A) signal.....	1213-1338
SV40 origin/promoter.....	1339-1903
Puromycin resistant gene	1904-2503
SV40 poly(A) signal	2665-2698
beta-lactamase (Amp ^R) gene	2807-3667
pUC replication origin (pUC Ori).....	3846-4650
Transcriptional terminator (Ter).....	4651-5258
Lac operator 1 (Lac O1)	5259-5279
Transcriptional pause site (TPS).....	5412-5483