


Format: 96 tip  
 Tip Volume Range Selection: 10 µL - 1000 µL  
 Material: Conductive LLS Black Polypropylene  
 Certified RNase-/DNase-free, nonpyrogenic  
 Available Options:  
 • With or without aerosol-resistant filters  
 • Nonsterile or sterilized



Tip Box Format: Non-ANSI hanging tip rack  
 Material: Number 7 Polycarbonate

**Axygen<sup>®</sup> Robotic Tips**

Type	Max Volume	Compatible Head	Tip Box	Packaging	Standard Cat. No.
 Conductive, Non-Filtered	20 µL	8-Channel		Nonsterile Sterile Quantity: 96 tips/rack, 24 racks/case	TT-20-CBK-HTR TT-20-CBK-HTR-S
 Conductive, Filtered	10 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-20-CBK-HTR-S
 Conductive, Non-Filtered	50 µL	8-Channel		Nonsterile Sterile Quantity: 96 tips/rack, 24 racks/case	TT-50-CBK-HTR TT-50-CBK-HTR-S
 Conductive, Filtered	25 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-50-CBK-HTR-S
 Conductive, Non-Filtered	200 µL	8-Channel		Nonsterile Sterile Quantity: 96 tips/rack, 24 racks/case	TT-200-CBK-HTR TT-200-CBK-HTR-S
 Conductive, Filtered	175 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-200-CBK-HTR-S
 Conductive, Non-Filtered	1000 µL	8-Channel		Nonsterile Sterile Quantity: 96 tips/rack, 16 racks/case	TT-1000-CBK-HTR TT-1000-CBK-HTR-S
 Conductive, Filtered	1000 µL	8-Channel		Sterile Quantity: 96 tips/rack, 16 racks/case	TTF-1000-CBK-HTR-S