

-25°C Freezer



400L

Technical Parameters

Model	BDF-25V400	
External Size(W*D*H) mm	600*645*1875	
Capacity	400L	
Type	Vertical	
Temp. Range	-10°C~-25°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor control	
Alarm	Sensor failure alarm; Open-door alarm; High or low temperature alarm; Power failure alarm and low battery alarm	
Refrigeration Type	Direct refrigeration	
Refrigerant	R134a, CFC free	
Condenser&Evaporator	Wire tube condenser	
Construction	Structure	Unibody design and thickened environmental protection foam insulation layer
	Internal Material	ABS Material
	External Material	Cold-rolled steel coated with anti-bacteria powder
Shelf/Drawer	6 wire shelves	
Consumption	220W	
Power Supply	AC220V±10%;50/60HZ;	
Package Size(W*D*H) mm	650*720*2015	
Gross Weight	89kg	



268L



328L



270L/350L

Technical Parameters

Model	BDF-25V268	BDF-25V328	BDF-25V270	BDF-25V350
External Size(W*D*H) mm	673*676*1630	673*676*1886	690*744*1631	690*744*1885
Capacity	268L	328L	270L	350L
Type	Vertical			
Temp. Range	-10°C~-25°C			
Temp. Accuracy	0.1°C			
Control System	Microprocessor Control & LED Display			
Alarm	High and low temperature; System & sensor failure; Door ajar; Power failure alarm; Thermostat alarm		High or low temperature alarm, Sensor failure alarm, power failure alarm, abnormal door opening and closing alarm	
Refrigeration Type	Direct Refrigeration			
Refrigerant	R290, CFC free			
Compressor	International famous compressor			
Condenser & Evaporator	High efficiency air-cooled condenser, made of copper		Bundy tube condenser and inner coil tube evaporator	
Structure	Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)			
Construction	Internal Material	Cold-rolled steel coated with anti-bacteria powder		SUS304
	External Material			Cold-rolled steel coated with bacterial powder
Shelve	3pcs	4pcs	3pcs (without inner door)	4pcs (without inner door)
Consumption	280W	280W	290W	305W
Power Supply	AC220V±10%, 50/60HZ		110V60HZ	220V±10%;50Hz
Package Size(W*D*H) mm	760*770*1800	770*770*2040	810*860*1820	810*860*2100
Gross Weight	125kg	140kg	155kg	170kg