

## Laboratory Refrigerator



656L



1500L

### Technical Parameters

Model	BPR-5V650	BPR-5V1500
External Size(W*D*H) mm	1220*642*1885	1800*748*1965
Capacity	656L	1500L
Temp. Range	2~8°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor Control & Large Screen LED Display	
Alarm	Audio and visual alarm for: High and low temperature, Thermostat fault, power failure(only for BPR-5V650), Door ajar	
Refrigeration Type	Forced air refrigeration system	
Refrigerant	R600a, CFC Free	R134a, CFC Free
Condenser&Evaporator	High efficiency air-cooled condenser, made of copper	
Defrost	Auto defrost	
Construction	Structure	Unibody design & Mono-assembly foaming (rigid polyurethane insulation material)
	Internal Material	Cold-rolled steel coated with anti-bacteria powder SS304
	External Material	Cold-rolled steel coated with anti-bacteria powder SS304
Shelves	Adjustable shelves made of high quality steel wire	
	10 pcs	12 pcs
Consumption	310W	830W
Power Supply	AC110/220V±10%, 50/60Hz	
Standard Accessory	LED lamp, Shelves, Door keys, Testing hole, Battery	
Package Size(W*D*H) mm	1290*695*1940	1940*880*2035
Gross Weight	188kg	190kg

## -25°C Freezer



BDF-25H305

### Features

1. Digital temperature display.
2. The safety lock can prevent random opening.
3. Baskets made of steel wire are convenient to store articles.
4. Round corner design of doors and the body provides luxury and elegance.
5. Delayed start and safe stop interval between restart and being terminated, ensure reliable running.

### Technical Parameters

Model	BDF-25H305	
External Size(W*D*H) mm	1000*520*665	
Capacity	305L	
Temp. Range	-10°C~-25°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor control&LED display	
Alarm	High or low temperature alarm, Temperature sensor failure alarm, over range alarm.	
Refrigeration Type	Direct refrigeration	
Refrigerant	R600a, CFC free	
Condenser&Evaporator	Internal-mounted condenser and coil type evaporator	
Construction	Structure	Whole foaming, freon cyclopentane/isopentane high density polyurethane foam
	Internal Material	Aluminum
	External Material	Cold-rolled steel coated with anti-bacteria powder
Basket	none	
Consumption	106W	
Power Supply	AC220V±10%,50HZ	
Package Size(W*D*H)mm	1223*777*1005	
Gross Weight	97kg	