


Format: 96 tip
 Tip Volume Range Selection: 20 µL - 1000 µL
 Material: Conductive LLS Black Polypropylene
 Certified RNase-/DNase-free, nonpyrogenic
 Available Options:
 • With or without aerosol-resistant filters
 • Nonsterile or sterilized

Tip Box Format: ANSI
 Material: Number 5 Polypropylene

1 mL Tip Box Format: Non-ANSI hanging tip rack
 Material: Number 7 Polycarbonate



Axygen® Robotic Tips

Type	Max Volume	Compatible Head	Tip Box	Packaging	Standard Cat. No.
Non-Filtered	20 µL	Variospan-1, 4, 8		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/box, 5 packs/case	RRP-20-CBK-R RRP-20-CBK-R-S
Non-Filtered	50 µL	Variospan-1, 4, 8		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/box, 5 packs/case	RRP-50-CBK-R RRP-50-CBK-R-S
Filtered	25 µL	Variospan-1, 4, 8		Sterile Quantity: 96 tips/rack, 10 racks/box, 5 packs/case	RRF-25-CBK-R-S
Non-Filtered	235 µL	Variospan-1, 4, 8		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/box, 5 packs/case	RRP-200-CBK-R RRP-200-CBK-R-S
Filtered	175 µL	Variospan-1, 4, 8		Sterile Quantity: 96 tips/rack, 10 racks/box, 5 packs/case	RRF-175-CBK-R-S
Non-Filtered	1000 µL	Variospan-1, 4, 8		Nonsterile Sterile Quantity: 96 tips/rack, 16 racks/case	TT-1000-CBK-HTR TT-1000-CBK-HTR-S
Filtered	1000 µL	Variospan-1, 4, 8		Sterile Quantity: 96 tips/rack, 16 racks/case	TTF-1000-CBK-HTR-S

