

Celleste™ System Requirements

Pub. No. MAN0025839 Rev. A.0

	Recommended	High Performance for 3D Module
Operating system	Windows™ 10 (64 bit) - Anniversary Update, Professional, and up	Windows™ 10 (64 bit) - Anniversary Update, Professional and up
	Windows™ 8.1 (64 bit) - Professional and up	
	Windows™ 7 (64 bit) - SP1, Professional, and up	
	32 bit versions are supported, but not recommended.	
Free disk space	Multi high-speed SATA hard disks or SSDs	Multi high-speed SATA hard disks or SSDs configured in RAID 5
	8 GB free on installation drive + free space for images (20 GB+)	16 GB free on installation drive + free space for images (1 TB+)
Processor	2.8 GHz CPU Intel™ quad-core processor or better	3.0 GHz CPU Intel™ quad-core processor or better
RAM	16 GB memory or higher	48 GB memory or higher
Graphics card	NVIDIA™ GeForce GTX™ cards with 4 GB graphics memory & OpenGL™ 4.2 or higher	NVIDIA™ GeForce GTX™ cards with 8 GB graphics memory & OpenGL™ 4.2 or higher
Drives	DVD-ROM drive (if optional DVD is purchased)	DVD-ROM drive (if optional DVD is purchased)
USB port	Required for hard licenses and online license activation	Required for hard licenses and online license activation
Internet connection	Internet connection required for online services ^[1]	Internet connection required for online services ^[1]
Internet browser	Internet Explorer™ version 9 or higher	Internet Explorer™ version 9 or higher

^[1] Online services: License activation/registration and product feedback/updates.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

thermofisher.com/support | thermofisher.com/askaquestion



Life Technologies Corporation | 29851 Willow Creek | Eugene, Oregon 97402 USA

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0025839

Revision	Date	Description
A.0	07 October 2021	New user bulletin for <i>Celleste™ System Requirements</i> .

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2021 Thermo Fisher Scientific Inc. All rights reserved. NVIDIA, NVIDIA GeForce, and NVIDIA GeForce GTX are trademarks and/or registered trademarks of NVIDIA Corporation. Intel is a trademark of Intel Corporation. Internet Explorer and Windows are trademarks of the Microsoft Corp. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For Research Use Only. Not for use in diagnostic procedures.