

Thermo Scientific Nalgene Centrifugeware Adapters Centrifuge Bottle Adapter

Cat. No. DS3125-0250

The adapter is designed for use with 250mL bottles in Thermo Scientific™ Sorvall™ GSA rotors used in Sorvall high-speed centrifuges. It can replace Sorvall adapter numbers 00424 and 00516. Unlike rubber adapters, the Thermo Scientific Nalgene™ adapter will not cling to bottles after spinning, and it makes bottles easier to remove from the rotor.

The Nalgene adapter fits most 250 mL centrifuge bottles, including Thermo Scientific Cat. Nos. 3120-0250, 3121-0250, 3122-0250, 3127-0250, 3140-0250 and 3141-0250.

Support for Spherical-Bottom Centrifuge Bottles
Cat. No. DS3124-0010 (white polycarbonate, autoclavable)

This support is designed for use with Nalgene 250 mL Spherical-Bottom Centrifuge Bottles, Cat. No. DS3123-0250.

It can be used with the following rotors: Sorvall H4000, H6000A, HG-4L, GSA; Beckman JS-3.0, JS-4.0, JS-4.2, JS-5.2, JS-7.5.

For optimal performance

Clean these adapters in a non-abrasive, non-alkaline cleanser and rinse thoroughly. Maintain the proper maintenance schedule for your rotors. Do not allow deposits to develop under the adapter.

For technical assistance, contact Technical Service at: **technical.nalgene@thermofisher.com**. Visit **www.thermoscientific.com** product resource pages for complete product use guidelines.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

8-0403-49 1212

Thermo
SCIENTIFIC

Contact us for Sales and Service
thermoscientific.com/contactus

*Contact information contained within
this document may be incorrect.