

**Catalog Number:** 431489

**Product Description:** Corning® Spin-X® UF Concentrator, 20 mL, 30k MWCO,  
Non-Sterile

**Component Materials:**

- Tube - Virgin Polycarbonate, meets *USP, Class VI* requirements for plastic containers and closures.
- Cap - Virgin High Density Polyethylene, Heavy metal free (meets *CONEG* req.) color concentrate.
- Membrane - Polyethersulfone (PES).

**Product Dimensions:**

- Tube Height w/Cap - 4.57 inches                      Cap OD - 1.38 inches
- Tube Diameter at Base - 1.18 inches

**Maximum Fill Volume and RCF:**

	<u>Max. Volume</u>	<u>Max. RCF</u>
Swing out (bucket)	20 mL	4,000xg
Fixed Angle (25°)	14 mL	6,000xg

**Membrane Retention:**

The product has been tested for retention using a standard molecular weight marker solution.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

- Integrity Test: Pressure hold test method.
- Visual Attributes: Visual examination of product.
- Packaging: Visual examination of accurate labeling and correct product configuration.

**Revision Date:**

03-14-19

Rev No: 3