

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 1300P-500

Product Description: Corning® Bottle, Washing, Reusable Plastic, LDPE, Not Sterile, 500 mL
GL-25 PP Cap and PP Spray Tube

Component Materials:

Bottle - Low Density Polyethylene (LDPE)
Stopper - Polypropylene (PP)
Spray Tube - Polypropylene (PP)

Attributes:

Capacity - 500 mL
Product Dimensions - 74 mm Body OD x 180 mm High without Cap

Autoclaving:

This product IS NOT suitable for Autoclaving.

Heating / Cooling:

The upper temperature limit for this plastic material is 80°C.

This product IS NOT suitable for Dry Heat Sterilization as it will melt. This product MUST NOT be heated by any direct, surface contacting heat source.

The lower limit usage temperature for this product is -50°C.

Vacuum / Pressure:

This product is not designed for and is not rated for vacuum or pressure applications.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

01 OCT 2012

Rev No: 1