

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 4519-120

Product Description: Media Handling Fitting, Vertical Sidearm Flask, Dual 1/8" - 1/4", 8 & 15 L

Component Materials:

Cap - Polypropylene (PP)
Nuts - Polyphenylene Oxide (PPO) / Polystyrene (PS)
Insert - Polyethylene terephthalate (PET)
O-ring - Synthetic Rubber
Tube - Stainless Steel

Attributes:

Product Dimensions - 54 mm OD (Cap) x 457mm High (Tube)

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

Heating / Cooling:

This product is primarily designed for a specific use that does not involve direct heating. In normal use, temperatures will be from -50° to +125°C. This is below the temperature at which thermal failure could be a problem.

Vacuum / Pressure:

This product is not designed for and is not rated for vacuum or pressure applications.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

5 AUG 2013

Rev No: 1