

## Forced Air Drying Oven

### Features:

- \* External material: cold-rolled steel with anti-bacteria powder coating.
- \* Stainless steel inner chamber, round angle structure, adjustable shelves.
- \* Double-layer glass observation window, adjustable airtight buckle lock.
- \* PID control with LED display.
- \* Over-temperature protection.



### Technical Parameters:

Model	BOV-T25F	BOV-T50F	BOV-T105F	BOV-T200F
Capacity	25L	54L	108L	212L
Temp. Range	50~200°C			
Temp. Precision	0.1°C			
Temp. Fluctuation	±1°C			
Ambient Temp.	5~40°C			
Timing Range	1~9999min			
Shelves No.	2 pcs			
Consumption	650W	1000W	1600W	2500W
Power Supply	Standard: 220V±10%, 50/60Hz; Optional: 110V, 60Hz			
Internal Size(W*D*H) mm	298*290*274	419*332*348	550*348*550	602*498*602
External Size(W*D*H) mm	580*554*455	708*595*532	830*620*735	880*775*785
Package Size(W*D*H) mm	710*660*620	840*690*700	960*720*900	1010*880*940
Gross Weight(kg)	51	66	115	145