

UV Lamp

Handheld UV Lamp

Application:

The handheld UV Lamp irradiate the single-cell microorganisms such as cells and viruses with UV light to destroy the structure of DNA in the center of their lives, so that the proteins that make up this organism cannot be synthesized, causing them to immediately die or lose their ability to reproduce. Under normal circumstances, UV light irradiation for more than 10 seconds can achieve the effect of sterilization. It is suitable for the sterilization and disinfection of computer keyboards, tablet PC, mobile phones, desktops, door handles and kitchen utensils.



BK-UVC-01A



BK-UVC-01B

Technical Parameters:

Model	BK-UVC-01A	BK-UVC-01B
Consumption	3W	1W
Light Tube Size	100mm	
Lamp Wavelength	253.7nm (no ozone)	
UV Lamp Service Life	6000 hours	
Material	Shell: ABS; Light tube: Quartz	
Operating Temperature	0~60°C	
Relative Humidity	10~75% RH	
Power Supply	4 AAA batteries, or connect to USB power	2 AAA batteries
External Size (H*W*D)	125*35*25mm	190*30*25mm
Net Weight	60g (No batteries)	50g (No batteries)
Inner Box Size (W*D*H)	156*96*30mm	190*30*25mm
Package Details	100 pcs / carton	
Package Size (W*D*H)	156*96*30mm (1 box)	190*30*25mm (1 box)
	495*325*320mm (1 carton)	345*285*210mm (1 carton)
Gross Weight	0.16kg (1 box)	0.05kg (1 box)
	17kg (1 carton)	6kg (1 carton)