Corning Incorporated Life Sciences

Registered ISO 9001:2008

Product Description

Catalog Number: 3680

Product Description: Corning 384-well, clear plate, with lid

Component Materials:

Plate/Lid - Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and

closures.

Product Dimensions:

Length of Plate - 5.030 in. Width of Square Well @ - .143 in.

Top

Width of Plate - 3.365 in. Width of Square Well @ - .105 in.

Bottom

Depth of Well - .450 in. Height with Lid - .637 in.

Tolerances of - +/-.010 in. Working Volume per well - 20 μL to 75 μL

Dimensions

Sterilization:

This lot has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization.*Sterility Assurance Level: SAL 10⁻³

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for seal and barrier integrity, accurate labeling, and correct

product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 5