



Format: 96 tip and 384 tip  
 Tip Volume Range Selection: 20 µL - 200 µL  
 Material: Clear Polypropylene  
 Certified RNase-/DNase-free, nonpyrogenic  
 Available Options:  
 • With or without aerosol-resistant filters  
 • Nonsterile or sterilized  
 • Standard or MAXIMUM RECOVERY®



Tip Box Format: ANSI  
 Material: Number 5 Polypropylene



**Axygen® Robotic Tips**

Type	Max Volume	Compatible Head	Tip Box	Packaging	Standard Cat. No.	MAXIMUM RECOVERY® Cat. No.
Non-Filtered 	20 µL	MDT-96 P50		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PK-20-R PK-20-R-S	PK-20-L-R PK-20-L-R-S
Filtered 	20 µL	MDT-96 P50		Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PKF-20-R-S	PKF-20-L-R-S
Non-Filtered 	50 µL	MDT-96 P50		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PK-50-R PK-50-R-S	PK-50-L-R PK-50-L-R-S
Filtered 	30 µL	MDT-96 P50		Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PKF-30-R-S	PKF-30-L-R-S
Non-Filtered 	235 µL	MDT-96 P235		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PK-200-R PK-200-R-S	PK-200-L-R PK-200-L-R-S
Filtered 	180 µL	MDT-96 P235		Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PKF-180-R-S	PKF-180-L-R-S
Wide-bore, Non-Filtered 	235 µL	MDT-96 P235		Nonsterile Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PK-205-WB-R PK-205-WB-R-S	
Wide-bore, Filtered 	180 µL	MDT-96 P235		Sterile Quantity: 96 tips/rack, 10 racks/pack, 5 packs/case	PKF-205-WB-R-S	
Non-Filtered 	30 µL	MDT 384-30 µL		Nonsterile Sterile Quantity: 384 tips/rack, 10 racks/pack, 5 packs/case	PK-384-R PK-384-R-S	
Filtered 	30 µL	MDT 384-30 µL		Sterile Quantity: 384 tips/rack, 10 racks/pack, 5 packs/case	PKF-384-R-S	