

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Product Description

Catalog Number: 9998-15

Product Description: PYREX® Screw Cap, Phenolic, PTFE Liner, G.P.I. 415-15

Component Materials:

Cap - Phenolic (Mineral Fill)

Liner - Polytetrafluoroethylene (PTFE) Faced Rubber

Attributes:

Product Dimensions - 18 mm OD x 16 mm High

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes. Caps should be loose on the container to avoid liners being pulled out by the vacuum created as the assembly cools.

Heating / Cooling:

This product is primarily designed for a specific use that does not involve direct heating. In normal use, temperatures will be from -50° to +125°C. This is below the temperature at which thermal failure could be a problem.

Vacuum / Pressure:

This product is a component part and not specifically rated for vacuum or pressure.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

16 AUG 2013

Rev No: 1