

PRP Centrifuge



Chamber

Features:

- * Max. speed 4500rpm.
- * 3 preloaded programs for fast operation.
- * 9 programs for customization can be stored in memory.
- * Powerful CPU for easy operation of speed, time and all other parameters.
- * Easy-to-read processing display and sound alert.
- * Automatic lid-lock release after running.
- * Maintenance-free brushless DC motor.

Technical Parameters:

Model	DM0412P
Speed Range	300-4500rpm, increment:10rpm
Max. RCF	2490×g
Timing Range	30sec-99min
Rotor Capacity	8*15ml, 12*10ml/7ml/5ml
Display	LCD
Preloaded Program	PRP, PRF, CGF
Custom Programs Memory	9 programs (6 steps for each program)
Noise	≤56dB
Power	Single-phase, 100V-240V,50Hz/60Hz, 3A,70W
External /Package Size (W*D*H)	301*354*217mm/460*420*320mm
Net/Gross Weight	6/9kg