

Catalog Number: 430499

Product Description: Corning ® cap inserts for cryogenic vials, assorted colors

Insert Colors:

Blue, green, yellow, red, white.

Component Materials:

Cap Insert - Virgin Polypropylene, meets *USP, Class VI* requirements for plastic containers and closures. Heavy metal free (meets CONEG req.) color concentrate.

Product Dimensions:

Diameter insert - 0.38 in.

Length of insert - 0.20 in.

Tolerances - +/- 0.05 in.

BSE/TSE:

This product complies with the latest revision of EMEA/410/01 "Note for Guidance on minimizing the risk of transmitting animal spongiform encephalopathy agents via human veterinary medicine products" by virtue of all bovine derived material having been processed per specific conditions of section 6.4 of EMEA/410/01.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Packaging: Inspection for accurate labeling and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 4