CORNING	Life Sciences	Document Number: LSR00137	Rev.: 5
Description (Class): Reference, Customer 1 35x10 mm Easy Grip D		Page: <b>1 of 2</b>	

# **APPLICATIONS:**

- Growth of cells.
- Unimpaired observation of dish contents.
- Ideal for microphotography.

### FEATURES:

- Flat, distortion free optics.
- Lids designed for optimal gas exchange.
- Stacking rings allow for easier stacking and handling.
- Frosted rim improves the handling of the dish. This feature facilitates the ability to pick up the dish without accidentally removing the lid to guarantee aseptic manipulation.
- Reliable vacuum-gas plasma treatment cleans and permanently and consistently modifies the cellgrowth surface.
- Non-pyrogenic
- Packaged in peel-open medical-style bags.
- Durable construction for stable dish manipulation

# **REGULATORY COMPLIANCE:**

Falcon<sup>®</sup> products are manufactured under the current ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

## STERILITY:

Product is gamma irradiated

### **RESTRICTIONS:**

• Refer to "Thermoplastics Properties Chart," available at: https://www.corning.com/catalog/cls/documents/compatibility-guides/CLS-DL-AN-226.pdf


Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

CORNING	Life Sciences	Document Number: LSR00137	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet 35x10 mm Easy Grip Dish			Page: <b>2 of 2</b>

# All dimensions in millimeters (mm)

- Assembled Dish and Lid height: 12.8 mm Effective growth area = 9.6 cm<sup>2</sup>

- Working Volume = 2.5 3.0 mL Actual dimensions: 40.28 mm O.D. x 6.17 mm
- Actual growth area: 11.78 cm<sup>2</sup>



#### MATERIAL: CRYSTAL-GRADE POLYSTYRENE

CAT. NO.		DESCRIPTION	DISH AND LID KEY DIMENSIONS											
		DESCRIPTION	Α	В	С	D	E	F	G	Н	J	Κ	L	М
1	353001	35x 10 mm Easy Grip Dish. TC-treated. 20/sleeve, 500/case		35.61	35.00	37.90	39.14	11.05	9.56	1.01	40.28	38.78	6.17	.889
2	351008	35x 10 mm Easy Grip Dish Non-TC-treated. 20/sleeve, 200/case												
2		353801 35x 10 mm Easy Grip Dish. Corning® Primaria™-reated. 20/sleeve,												
		500/case		i				i						

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners. © 2013-2019 Corning Incorporated. All rights reserved.