 <b>Life Sciences</b>	<b>Affected Site(s): US : MA : CMO Tewksbury; US : MA : Tewksbury</b>	Document Number: <b>LSR00220</b>	
		Revision: <b>1</b>	Page: <b>1 of 2</b>
Description: <b>Falcon 3mL Transfer Pipet</b>			

**FALCON® 3mL Transfer Pipets**

**FEATURES:**

Two-styles of packaging- either bulk-style (non-sterile) or individually wrapped (sterile) are available.

Graduations at 1 & 2 mL with total capacity of 3mL

Provide quick and easy liquid transfers

**CLAIMS:**

Quality Management System - Complies with the current version of the ISO 9001 Standard.

USP Class VI Testing - All material resin is tested, qualified and shown to be non-toxic as established in the Standards USP Class VI Chapter, "Biological reactivity Tests, in Vitro" and Chapter, "Biological Reactivity Tests, in vivo".

Sterility (**#357575 only**) - Product has been sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137, "Sterilization of health care products- Radiation". Products meet a minimum Sterility Assurance Level (SAL) of 10<sup>-3</sup>.

Quality Control Testing - Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release. - This product met Corning Incorporated - Life Sciences' high standards of quality at the time of batch/lot release.

**PRODUCT AND PACKAGING CONFIGURATION**

Catalogue #	Capacity	Sterility	Case Qty	Pack Qty	Packaging Style
357575	3mL	SAL 10 <sup>-3</sup>	500 pcs	1/pk	Individually Wrapped
357524	3mL	Non-sterile	1000 pcs	500/pk	Bulk-style

**MATERIALS**


Low-density polyethylene (LDPE)

**LOT NUMBER DESIGNATION (7-DIGITS):**


First digit – Last number of year of manufacture; next 3 digits – Julian date, start of manufacturing; last 3 digits – batch identification

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks). All other trademarks are the property of their respective owners. © 2020 Corning Incorporated. All rights reserved.

 <b>Life Sciences</b>	Affected Site(s): <b>US : MA : CMO  Tewksbury; US : MA :  Tewksbury</b>	Document Number: <b>LSR00220</b>	
		Revision: <b>1</b>	Page: <b>2 of 2</b>
Description: <b>Falcon 3mL Transfer Pipet</b>			

**PRODUCT DIMENSIONS**



**\*All dimensions in millimeters**

**RESTRICTIONS:**

Refer to “Thermoplastics Properties Chart,” (CLS-DL-AN-226) available at [www.corning.com/lifesciences](http://www.corning.com/lifesciences):  
<https://www.corning.com/catalog/cls/documents/compatibility-guides/CLS-DL-AN-226.pdf>

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks). All other trademarks are the property of their respective owners. © 2020 Corning Incorporated. All rights reserved.