



The following features are present in the vector based on the nucleotide sequence.

Cytomegalovirus (CMV) promoter.....	8-635
<i>Gaussia-Dura</i> luciferase gene	646-1203
BGH poly(A) signal.....	1212-1337
SV40 origin/promoter.....	1338-1902
Puromycin resistant gene	1903-2502
SV40 poly(A) signal	2664-2697
beta-lactamase (Amp ^R) gene	2806-3666
pUC replication origin (pUC Ori).....	3845-4649
Transcriptional terminator (Ter).....	4650-5257
Lac operator 1 (Lac O1)	5258-5278
Transcriptional pause site (TPS).....	5411-5482