



The following features are present in the vector based on the nucleotide sequence.

Thymidine kinase (TK) promoter	26-777
<i>Gaussia-Dura</i> luciferase gene	790-1347
BGH poly(A) signal.....	1356-1481
SV40 origin/promoter.....	1482-2046
Puromycin resistant gene	2047-2646
SV40 poly(A) signal	2808-2841
beta-lactamase (Amp ^R) gene	2950-3810
pUC replication origin (pUC Ori).....	3989-4793
Transcriptional terminator (Ter)	4794-5401
Lac operator 1 (Lac O1)	5402-5422
Transcriptional pause site (TPS).....	5555-5626