Corning Incor Life Scien <u>Registere</u> ISO 9001:20	Product Description
Catalog Number:	3923
Product Description:	Costar® 1 x 8 Stripwell <sup>TM</sup> , high binding, solid white plate, flat bottom, without lid
Component Materials: Strips-Holder-	Virgin Polystyrene, meets <i>USP, Class VI</i> requirements for plastic containers and closures. White concentrate Polystyrene, Heavy metal free (meets <i>CONEG</i> req.) color concentrate.
<b>Product Dimensions:</b> Length of Plate Width of Plate Depth of Well Tolerances of Dimensions	-5.025 in.Diameter of Well @ top270 in3.370 in.Diameter of Well @250 in420 in.Height without Lid560 in+/010 in.Working Volume per well36 mL

## Surface Characterization:

Surface is characterized to be high binding.

## **Optical Characteristics:**

The product is made of opaque white polystyrene to minimize well to well crosstalk and background fluorescence and/or luminescence.

### Binding Characteristics:

Each lot has been tested for the attribute of well to well protein binding utilizing a luminescent immunoassay to measure precision. Coefficient of variation must be  $\leq 5\%$ .

### Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.Visual Attributes:Visual examination of the product.Packaging:Inspection for accurate labeling, and correct product configuration.Functional:Luminescent immunoassay.

# Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 5