

# BIOLGX<sup>®</sup> NEW Instruments

## Liquid Handling, Centrifuges, Mixers & Thermocyclers

### Mixers & Thermocyclers



01-120X (MX-RD-E)



01-130X (MX-RL-E)

01-620X (SK-O180-E)  
01-610X (SK-L180-E)  
01-630X (SK-R1807-E)  
01-640X (SK-D1807-E)



01-700X (TC1000-G)



01-110X (MX-S)  
With standard top cup VT1.1 adjustable speed

### Magnetic Stirrers



01-320X  
Ultra-flat Compact Magnetic Stirrer



01-310X (MS-H280-Pro)  
LED Digital Magnetic Hotplate Stirrer

### Palm Micro Centrifuges

Small Machines, Multiple Applications



CE FC

01-800X (D1008)  
CE explosion-proof safe test  
IEC/EN6 1010-2-20

### Liquid Handling

Comfortable  
Finger support

Simply turn the plunger  
button for volume selection

Tip ejector allows  
convenient one-handed  
operation

Ejector collar and tip  
cone can be removed

Durable tip cone  
materials provide  
excellent chemical  
resistance



**MicroPette plus  
8-channel Pipettes**  
01-2801 ( 0.5-10 $\mu$ L )  
01-2805 ( 5-50 $\mu$ L )  
01-2830 ( 50-300 $\mu$ L )

#### Plus Autoclavable Pipettes 01-220X

- Easy for one-handed operation
- Pipette controller for use with pipettes from 0.1 - 100mL
- Powerful motor fills a 25mL pipette in <5 seconds
- Supplied with integral but replaceable 0.45 $\mu$ l filter

The dispense button is pressed only up to a point where a check can be felt, the dispensing function will be based on the force of gravity and the liquid will simply flow out.



Simply turn the plunger  
button for volume selection

Tip ejector allows  
convenient one-  
handed operation

**MicroPette plus  
Single-channel Pipettes**  
01-2101 ( 0.5-10 $\mu$ L )  
01-2102 ( 20-200 $\mu$ L )  
01-2110 ( 100-1000 $\mu$ L )