

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 2947-75X38

Product Description: Corning® Microscope Slide, Plain, 75 x 38 mm

Component Materials:

Type II - Soda-Lime Glass, per ASTM E-438 (Standard Specification for Glasses in Laboratory Apparatus)

Attributes:

Product Dimensions - 75 mm x 38 mm x 0.90 - 1.10 mm Thick

Autoclaving:

This product is made of Soda Lime glass and IS NOT designed to be autoclaved.

Heating / Cooling:

This product is designed for use in general laboratory microscopy. It is not intended for heating or cooling.

Vacuum / Pressure:

This product is not designed for and is not rated for vacuum or pressure applications.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1