

## Mould Incubator



LCD Display



USB Port

### Application:

This series of products are widely used in medical and health, epidemic prevention and inspection, biopharmaceuticals, animal husbandry research, environmental protection, agricultural scientific research and other fields.

### Features:

- \* LCD display.
- \* The refrigeration adopts a new fluorine-free design to ensure green and environmental protection.
- \* With parameter memory function, the call will automatically resume operation.
- \* The inner box is designed with a circular arc structure, which is easy to clean.
- \* Pullable movable shelves with adjustable spacing.
- \* The inner box is equipped with a fan to form a micro air circulation to improve the temperature uniformity in the box.
- \* Outer box material: cold-rolled steel plate sprayed on the surface; internal material: SUS304 stainless steel plate.

### Technical Parameters:

Model	BJPX-M80N	BJPX-M100N	BJPX-M150N	BJPX-M200N	BJPX-M250N	BJPX-M300N	BJPX-M400N
Capacity	80	100	150	200	250	300	400
Temp.Control Range	0°C~65°C(Microcomputer PID control)						
Humidity Range	50%~90%RH(Microcomputer PID control)						
Temp.Uniformity	±2°C						
Temp.Fluctuation	±1°C						
Humidity Fluctuation	±5~8%RH						
Temperature and Humidity Display Accuracy	0.1°C / 1%RH						
Timing Range	0~9999min(h)						
Power Consumption	1200W	1200W	1200W	1300W	1300W	1500W	1500W
Refrigerant	R134a					R404a	
Power	AC220V, 50/60Hz(Standard); AC110V,60Hz(Optional)						
Accessories	Standard: Socket, lighting; Optional:USB interface, Printer						
External Size(W*D*H)	540*580 *1110mm	540*580 *1230mm	590*630 *1390mm	590*660 *1540mm	620*700 *1610mm	660*700 *1710mm	720*760 *1810mm
Internal Size(W*D*H)	400*370 *560mm	400*370 *680mm	450*420 *850mm	450*450 *1000mm	480*490 *1070mm	520*500 *1170mm	580*540 *1270mm
Package Size(W*D*H)	670*710 *1280mm	670*710 *1400mm	720*760 *1560mm	720*790 *1710mm	750*830 *1780mm	790*830 *1880mm	850*880 *1980mm
Net/Gross Weight	55/90kg	60/95kg	100/130kg	105/148kg	113/157kg	125/167kg	145/180kg

\* Notice:It cannot be used for long-term refrigeration and storage.