

(Automatic) Burettes

Standard type, compact & micro burettes

Capacity [ml]	Subdivision [ml]	Height [mm]	Inner diameter [mm]	Error limit [± ml]	
				BLAUBRAND®	SILBERBRAND
2	0.01	600/660*	3.4	0.01	–
5	0.01	800/900*	3.4	0.01	–
10	0.02	800/900*	5	0.02	0.03
10	0.05	470	8	–	0.05
25	0.05	800	8	0.03	0.05
25	0.1	550	11.2	–	0.08
50	0.1	800	11.2	0.05	0.10

* Greater height only for micro burettes, Bang pattern, lateral stopcock incl. base

Automatic Burettes Dr. Schilling Pattern incl. bottle

Capacity [ml]	Subdivision [ml]	Height [mm]	Inner diameter [mm]	Error limit [± ml]
				SILBERBRAND
10	0.05	590	8	0.05
15	0.1	500	11.2	0.08
25	0.1	675	11.2	0.08
50	0.1	930	11.2	0.10

Automatic Burettes Pellet Pattern incl. bottle

Capacity [ml]	Subdivision [ml]	Height [mm]	Inner diameter [mm]	Error limit [± ml]	
				BLAUBRAND®	SILBERBRAND
10	0.02	1000	5	0.02	0.03
25	0.05	1000	8	0.03	0.05
25	0.1	700	11.2	–	0.08
50	0.1	1000	11.2	0.05	0.10

Compact Automatic Burettes

Capacity [ml]	Subdivision [ml]	Height [mm]	Inner diameter [mm]	Error limit [± ml]	
				BLAUBRAND®	SILBERBRAND
10	0.02	790	5	0.02	–
10	0.05	470	8	–	0.05
25	0.05	790	8	0.03	–
25	0.1	550	11.2	–	0.08
50	0.1	790	11.2	0.05	0.10

