

TECHNICAL SHEET

PRODUCT NAME: Anti-Enterovirus pAb

Code: MT-20EV

Clone: NA

Specie: Rabbit pAb

Physical State: Liquid

Host Animal (source): Rabbit pAb

Specificity: 20EV is specific for *Enterovirus* detection.

Cross Reactivity: N/A

Purification: Anti-Strep A pAb has been purified in a multi-step procedure.

≥95% pure by SDS-PAGE.

Storage Conditions: Storage at 2-8°C.

Health & Safety Information: This product contains sodium azide as preservative. Although the

amount is lower than 0.09%, this substance is poisonous and

hazardous and it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance

with the respective national regulations.

Presentation: PBS pH: 7.4, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 18/01/2019



www.certest.es