

Thermo Scientific Nalgene Unwire and Microtube Racks

Thermo Scientific Nalgene ResMer manufacturing technology adds improved features and benefits to both Unwire™ and Microtube rack products. ResMer™ manufactured racks have excellent mechanical strength, and environmentally friendly formulation, and better heat and chemical resistance. Plus, colors stay bright longer, even after repeated processing.

Sterilization

ResMer manufactured racks can be sterilized by autoclaving for 20 minutes at 121°C, or by gas, radiation, or chemical disinfectants.

©2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

www.nalgenelabware.com

Asia: China Toll-free: 800-810-5118 or 400-650-5118; India: +91 22 6716 2200,
India Toll-free: 1 800 22 8374; Japan: +81 3 3816 3355 Other Asian countries: 65 68729717
Europe: Austria: +43 1 801 40 0; Belgium: +32 2 482 30 30; Denmark: +45 4631 2000;
France: +33 2 2803 2180; Germany: +49 6184 90 6940, Germany Toll-free: 08001-536 376;
Italy: +39 02 02 95059 or 434-254-375; Netherlands: +31 76 571 4440; Nordic/Baltic countries: +358 9 329 100

ResMer Chemical Resistance Summary	
Acids	Good
Bases	Good
Glycols	Good
Oils	Good
Chlorinated Hydrocarbons	Poor
Aromatic Hydrocarbons	Fair
Aliphatic Hydrocarbons	Good
Gasoline	Poor
Alcohols	Good
Salt Solutions	Good
Water	Good

Europe: Russia/CIS: +7 (812) 703 42 15; Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 12;
UK/Ireland: +44 870 609 9203
North America: USA/Canada +1 585 586 8800; USA Toll-free: 800 625 4327
South America: USA sales support: +1 585 899 7198
Countries not listed: +49 6184 90 6940 or +33 2 2803 2180

8-0407-93 0609

NALGENE®

Thermo
SCIENTIFIC